

Excellence In X-Ray Imaging

www.alerio.in

# OURO

**ALERIO**<sup>®</sup>  
ULTRAPORTABLE X-RAY SYSTEM



*Ai* ENABLED

## ULTRAPORTABLE X-RAY SYSTEM

For Accurate, Safe & Swift On-site TB Diagnosis



Introducing the All New **ALERIO OURO**  
Ultraportable light weight handheld X-ray  
system with integrated DR & AI Technology.



intertek



### SALIENT FEATURES

- > Laser guided X-Ray Field
- > X-Ray Field Size 17" X 17" at 1m SID
- > Tripod provision & Yoke
- > Intuitive Control Panel

### SAFETY FEATURES

- > Low Leakage Radiation
- > Active Over-Load Protection
- > Digital Close Loop Control



Introducing the All New **ALERIO OURO** Ultraportable light weight handheld X-ray system with integrated AI Technology.

## **Ai** ENABLED ULTRAPORTABLE X-RAY SYSTEM

For Accurate, Safe & Swift On-site TB Diagnosis

### GENERAL SPECIFICATION

Output kV Range	40 - 90 kV (1kV Step)
Output mA Range	8mA & 10mA
Focal Spot	0.8mm IEC Stationary Anode
Exposure Time Range	0.02 to 2 Sec (0.01 sec step)
Generator Type	High Frequency > 100kHz
Supply Voltage Range	180-250 VAC, 5 Amps 50/60 Hz
Weight	3.5 Kg
Controls	Microcontroller -closed loop system
User Interface	Graphical 2.4" TFT Display
Exposure Switch	Wired & Wireless remote with 2 stage triggers

